

# Hearing Aid Redesign: Test Plans

## ELECTRICAL TESTING

### Table of Contents:

Number	Test	Page
EE 1	Switch	2
EE 2	Speaker	7
EE 3	Sound Processing	11
EE 4	Microphone	14
EE 5	Battery Charger	18
EE 6	Bandpass and Pre-Amplification	20
EE 7	System	22

### Specification Summary:

Spec	Specification	Dir.	Units	Marginal	Ideal	Measured Value	Pass/Fail
S1	Lowest level of sound detected	min	dB	<60	<30		
S2	Frequencies amplified	max	Hz	300-7000	85-8000		
S3	Maximum amplification	target	dB	85-95	90		
S4	Levels of volume adjustment	max	levels	5	20		
S5	Time to charge earpiece	min	minutes	<60	<30		
S7	Connects to standard USB 2.0	target	yes/no	yes/no	yes		
S12	Battery life at max. Amplification	max	hours	>16	>48		

# Hearing Aid Redesign: Test Plans

## EE 1: Switch Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal
S4	Levels of Volume Adjustment	20	5

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

- \_\_\_ Switch
- \_\_\_ MC56F8006 Evaluation Kit
- \_\_\_ USB power cable
- \_\_\_ USB TAP cable

### Sections

- Part 1 – Test Volume Up and Down Interrupt
- Part 2 – Test Next and Previous Mode Interrupt
- Part 3 – Test Standby Interrupt
- Part 4 – Test Functionality of each Button Press

# Hearing Aid Redesign: Test Plans

## EE 1: Switch Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Test Volume Up and Down

- \_\_\_ 1. Solder wires to the switch
- \_\_\_ 2. Connect the volume up and down switches to the designated input pins on the eval board
- \_\_\_ 3. Flash the switch test software onto the Eval Board
- \_\_\_ 4. Press the volume up button
- \_\_\_ 5. Observe the volume up LED turns on and off when the button is pressed
- \_\_\_ 6. Press the volume down button
- \_\_\_ 7. Observe the volume down LED turns on and off when the button is pressed

### Summary of Data

Did the volume switches correctly interact with the eval board? Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 1: Switch Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 2 – Test Mode Next and Previous

- \_\_\_ 1. Solder wires to the switch
- \_\_\_ 2. Connect the mode next and previous switches to the designated input pins on the eval board
- \_\_\_ 3. Flash the switch test software onto the Eval Board (If needed)
- \_\_\_ 4. Press the mode next button
- \_\_\_ 5. Observe the mode next LED turns on and off when the button is pressed
- \_\_\_ 6. Press the mode previous button
- \_\_\_ 7. Observe the mode previous LED turns on and off when the button is pressed

### Summary of Data

Did the mode switches correctly interact with the eval board? Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 2 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 1: Switch Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 3 – Test Standby

- \_\_\_ 1. Solder wires to the switch
- \_\_\_ 2. Connect the standby switch to the designated input pins on the eval board
- \_\_\_ 3. Flash the switch test software onto the Eval Board (If needed)
- \_\_\_ 4. Press the standby button
- \_\_\_ 5. Observe the standby LED turns on and off when the button is pressed

### Summary of Data

Did the standby switch correctly interact with the eval board? Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 3 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 1: Switch Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 3 – Test Switch Functionality

- \_\_\_ 1. Solder wires to the switch
- \_\_\_ 2. Connect all of the switches to the evaluation board
- \_\_\_ 3. Flash the hearing aid software to the evaluation board
- \_\_\_ 4. Press the volume up and down buttons
- \_\_\_ 5. While debugging the code observe that the volume variable goes up and down
- \_\_\_ 6. Press the mode next and previous
- \_\_\_ 7. While debugging the code observe the mode variable changes accordingly
- \_\_\_ 8. Press the standby button
- \_\_\_ 9. While debugging the code observe that the code appropriately goes into standby mode

### Summary of Data

Did the switch function as intended? Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 4 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 2: Speaker Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal
S2	Frequencies Amplified	85 – 8000 Hz	300 – 7000 Hz

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

- Speaker
- Breadboard
- Signal Generator
- Signal Generator Test Leads

### Sections

- Part 1 – Determine Min Input to Speaker
- Part 2 – Speaker Functionality: Lowest Frequency Heard
- Part 3 – Speaker Functionality: Highest Frequency Heard

# Hearing Aid Redesign: Test Plans

## EE 2: Speaker Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Determine Min Input to Speaker

- \_\_\_ 1. Solder wires to the speaker
- \_\_\_ 2. Connect the speaker to VDD and GND
- \_\_\_ 3. Connect the data line of the speaker to the signal generator
- \_\_\_ 4. Set signal generator to a mid range frequency 1-2 kHz and lowest amplitude
- \_\_\_ 5. Observe if the speaker is outputting a tone, or sound
- \_\_\_ 6. If no sound heard up the amplitude a slight amount,  $x = \underline{\hspace{2cm}}$ ,
- \_\_\_ 7. Keep repeating steps 5 and 6 until a sound can just be heard
- \_\_\_ 8. Vary frequency to observe the speaker output at this level but different frequencies

### Summary of Data

Frequency (Hz)	Amplitude ()

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_



# Hearing Aid Redesign: Test Plans

## EE 2: Speaker Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 2 – Speaker Functionality: Lowest Frequency Heard

- \_\_\_ 1. Solder wires to the speaker (Ensure this is correct from previous part)
- \_\_\_ 2. Connect the speaker to VDD and GND (Ensure this is correct from previous part)
- \_\_\_ 3. Connect the data line of the speaker to the Signal Generator
- \_\_\_ 4. Set signal generator to a mid range value 1-3 kHz and the amplitude found above
- \_\_\_ 5. Observe the speaker output
- \_\_\_ 6. Lower the frequency of the signal generator
- \_\_\_ 7. Repeat steps 5 and 6 until the speaker no longer produces an audible noise.

### Summary of Data

Lowest Frequency Output by Speaker \_\_\_\_\_

Testing Part 2 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 2: Speaker Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 3 – Speaker Functionality: Highest Frequency Heard

- \_\_\_ 1. Solder wires to the speaker (Ensure this is correct from previous part)
- \_\_\_ 2. Connect the speaker to VDD and GND (Ensure this is correct from previous part)
- \_\_\_ 3. Connect the data line of the speaker to the Signal Generator
- \_\_\_ 4. Set signal generator to a mid range value 1-3 kHz and the amplitude found above
- \_\_\_ 5. Observe the speaker output
- \_\_\_ 6. Raise the frequency of the signal generator slowly
- \_\_\_ 7. Repeat steps 5 and 6 until the maximum rating of the speaker is reached or the ideal frequency is reached.

### Summary of Data

Highest Frequency Output by Speaker \_\_\_\_\_

Testing Part 3 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 3: Sound Processing Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal
S2	Frequencies Amplified	85 – 8000 Hz	300 – 7000 Hz
S3	Maximum Amplification	90 dB	85-95 dB

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

- \_\_\_ MC56F8006 Evaluation Kit
- \_\_\_ USB power cable
- \_\_\_ USB TAP cable
- \_\_\_ Matlab

### Sections

- Part 1: Test sound processing algorithm on known signal from Matlab
- Part 2: Test sound processing with signal generator

# Hearing Aid Redesign: Test Plans

## EE 3: Sound Processing Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Test sound processing with signal from Matlab

- \_\_\_ 1. Hard code signal into software
- \_\_\_ 2. Flash code onto evaluation board
- \_\_\_ 3. Set a break point after the signal is finished being processed before being output
- \_\_\_ 4. Run the code
- \_\_\_ 5. Take and plot the output of the sound processing algorithm in Matlab
- \_\_\_ 6. Compare the result to the initial signal

### Summary of Data

Did the sound processing algorithm modify the signal? Yes \_\_\_\_\_ No \_\_\_\_\_

Was the signal Amplified? Yes \_\_\_\_\_ No \_\_\_\_\_

Are modifications to the algorithm needed? Yes \_\_\_\_\_ No \_\_\_\_\_

Place Wave forms here

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 3: Sound Processing Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 2 – Test Sound Processing with Signal Generator

- \_\_\_ 1. Attach signal generator to the ADC ports of the eval board
- \_\_\_ 2. Flash code onto evaluation board
- \_\_\_ 3. Set a break point after the signal is finished being processed before being output
- \_\_\_ 4. Run the code
- \_\_\_ 5. Observe that the ADC interrupt LED goes on when interrupt is generated.
- \_\_\_ 6. Take and plot the output of the sound processing algorithm in Matlab
- \_\_\_ 6. Compare the result to the initial signal

### Summary of Data

Did the sound processing algorithm modify the signal? Yes \_\_\_\_\_ No \_\_\_\_\_

Did the ADC interrupt Trigger Yes \_\_\_\_\_ No \_\_\_\_\_

Was the signal Amplified? Yes \_\_\_\_\_ No \_\_\_\_\_

Are modifications to the algorithm needed? Yes \_\_\_\_\_ No \_\_\_\_\_

Place Wave forms here

Testing Part 2 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 4: Microphone Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal
S1	Lowest Level of Sound Detected	< 30 dB	<60 dB

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

- Directional Microphone
- Omni-Directional Microphone
- Breadboard
- Oscilloscope
- Oscilloscope Probes

### Sections

- Part 1 – Microphone Functionality: Normal speaking voice
- Part 2 – Microphone Functionality: Whisper
- Part 3 – Measure the output of the microphone

# Hearing Aid Redesign: Test Plans

## EE 4: Microphone Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Microphone Functionality: Normal speaking voice Test

- \_\_\_ 1. Solder wires to the microphones
- \_\_\_ 2. Connect the microphones to VDD and GND
- \_\_\_ 3. Connect the data lines of the microphones to the Oscilloscope
- \_\_\_ 4. Talk at a normal speaking level (approx. 60 dB)
- \_\_\_ 5. Observe oscilloscope output to obtain a quantifiable measurement of the microphones
- \_\_\_ 6. Save capture of oscilloscope to verify that it can pick up normal speech

### Summary of Data

Place capture of microphone waveforms here to prove functionality at a normal speaking level

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 4: Microphone Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 2 – Microphone Functionality: Whisper Test

- \_\_\_ 1. Solder wires to the microphones (Ensure this is correct from previous part)
- \_\_\_ 2. Connect the microphones to VDD and GND (Ensure this is correct from previous part)
- \_\_\_ 3. Connect the data lines of the microphones to the Oscilloscope
- \_\_\_ 4. Talk at a whisper(**approx. 30 dB**)
- \_\_\_ 5. Observe oscilloscope output to obtain a quantifiable measurement of the microphones
- \_\_\_ 6. Save capture of oscilloscope to verify that it can pick up speech at a whisper

### Summary of Data

Place capture of microphone waveforms here to prove functionality at a whisper speaking level

Testing Part 2 Sign Off \_\_\_\_\_ Date \_\_\_\_\_



# Hearing Aid Redesign: Test Plans

## EE 4: Microphone Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 3 – Measure the output of the microphone

- \_\_\_ 1. Solder wires to the microphones (Ensure this is correct from previous part)
- \_\_\_ 2. Connect the microphones to VDD and GND (Ensure this is correct from previous part)
- \_\_\_ 3. Connect the data lines of the microphones to the Oscilloscope
- \_\_\_ 4. Talk at varying levels of speech
- \_\_\_ 5. Observe oscilloscope output to obtain a quantifiable measurement of the microphones
- \_\_\_ 6. Save capture of oscilloscope to verify that it can pick up speech at a whisper
- \_\_\_ 7. Record the max voltage levels of the microphones output to determine the operating range of the microphone. This is useful to not harm electronics later on in the system

### Summary of Data

Place capture of microphone waveforms here to prove functionality at a whisper speaking level

Maximum voltage reading: Normal Level: \_\_\_\_\_

Maximum voltage reading: Low Level \_\_\_\_\_

Maximum voltage reading: High Level \_\_\_\_\_

Testing Part 3 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 5: Battery Charger Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal
S5	Time to charge Earpiece	<30	<60
S7	Connects to Standard USB 2.0 Port	Yes	Yes/No

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

- \_\_\_ Battery Charger IC
- \_\_\_ Batteries
- \_\_\_ Breadboard
- \_\_\_ Multimeter

### Sections

- Part 1: Test Battery Charger for Correct Operation

# Hearing Aid Redesign: Test Plans

## EE 5: Battery Charger Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Test Battery Charger for Correct Operation

- \_\_\_ 1. Wire up the charger circuit on the breadboard
- \_\_\_ 2. Connect the batteries
- \_\_\_ 3. Double check that the circuit is correct
- \_\_\_ 4. Apply the correct input voltage to the circuit
- \_\_\_ 5. Probe the circuit with the multimeter checking test point voltages
- \_\_\_ 6. Measure the batteries during charging to see if they are in fact charging
- \_\_\_ 7. Discharge the batteries to simulate circuit operation.
- \_\_\_ 8. Repeat steps 4 – 7 for multiple charge cycles

### Summary of Data

Did the battery charging circuit charge the batteries? Yes \_\_\_\_\_ No \_\_\_\_\_

Were we able to optimize the circuit at all? Yes \_\_\_\_\_ No \_\_\_\_\_

If so what was optimized:

What was the time to charge the earpiece? \_\_\_\_\_

Are we able to charge over USB 2.0? Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 6: Bandpass and Pre-Amplification Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### Sections

- Part 1: Test the Bandpass Filter and Pre-Amplifier

# Hearing Aid Redesign: Test Plans

## EE 6: Bandpass and Pre-Amplification Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Test Battery Charger for Correct Operation

\_\_\_ 1.

\_\_\_ 2.

\_\_\_ 3.

\_\_\_ 4.

\_\_\_ 5.

\_\_\_ 6.

\_\_\_ 7.

\_\_\_ 8.

### Summary of Data

Did the battery charging circuit charge the batteries? Yes \_\_\_\_\_ No \_\_\_\_\_

Were we able to optimize the circuit at all? Yes \_\_\_\_\_ No \_\_\_\_\_

If so what was optimized:

What was the time to charge the earpiece? \_\_\_\_\_

Are we able to charge over USB 2.0? Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Specifications Tested

Specification	Description	Ideal	Marginal
S1	Lowest Level of Sound Detected	< 30	< 60
S2	Frequencies Amplified	85 – 8000 Hz	300 – 7000 Hz
S3	Maximum Amplification	90 dB	85-95 dB
S4	Levels of Volume Adjustment	20	5
S7	Connects to USB 2.0 Port	Yes	Yes/no
S5	Time to Charge Earpiece	< 30	< 60
S12	Earpiece Battery Life at Max Amplification	> 48	> 16

### Revision History

Revision	Description	Date
1	Document Created	03/14/13

### Equipment

- \_\_\_ MC56F8006 Evaluation Kit
- \_\_\_ USB power cable
- \_\_\_ USB TAP cable
- \_\_\_ Directional Microphone
- \_\_\_ Omni-Directional Microphone
- \_\_\_ Speaker
- \_\_\_ Battery Charger Circuit
- \_\_\_ Switch
- \_\_\_ Filtering and Amplification Circuit

### Sections

- Part 1: Test for S1
- Part 2: Test for S2
- Part 3: Test for S3
- Part 4: Test for S4
- Part 5: Test for S7
- Part 6: Test for S5
- Part 7: Test for S12

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 1 – Test for S1

- \_\_\_ 1. Flash the latest code to the evaluation board
- \_\_\_ 2. connect all of the subsystems together
- \_\_\_ 3. Ensure batteries are charged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 4. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 5. Place the microphones and speaker inside of the audiologist test equipment
- \_\_\_ 6. Run a low sound scenario on the audiologist test equipment
- \_\_\_ 7. Compare the result of the low sound scenario to that of a known good hearing aid

### Summary of Data

Were we able to detect sound in the low sound scenario? Yes \_\_\_\_\_ No \_\_\_\_\_

Are modifications to the algorithm needed? Yes \_\_\_\_\_ No \_\_\_\_\_

Place Result Graph he

Testing Part 1 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 2 – Test for S2

- \_\_\_ 1. Flash the latest code to the evaluation board
- \_\_\_ 2. connect all of the subsystems together
- \_\_\_ 3. Ensure batteries are charged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 4. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 5. Place the microphones and speaker inside of the audiologist test equipment
- \_\_\_ 6. Run a scenario that will test at the limits of the frequency ranges to see hearing aid response
- \_\_\_ 7. Compare the result of the scenario to that of a known good hearing aid

### Summary of Data

Were we able to amplify the necessary frequencies? Yes \_\_\_\_\_ No \_\_\_\_\_

Are modifications to the algorithm needed? Yes \_\_\_\_\_ No \_\_\_\_\_

Place Result Graph here

Testing Part 2 Sign Off \_\_\_\_\_ Date \_\_\_\_\_



# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 3 – Test for S3

- \_\_\_ 1. Flash the latest code to the evaluation board
- \_\_\_ 2. connect all of the subsystems together
- \_\_\_ 3. Ensure batteries are charged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 4. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 5. Place the microphones and speaker inside of the audiologist test equipment
- \_\_\_ 6. Set mode to max gain
- \_\_\_ 7. Run a scenario and observe the maximum gain of the hearing aid
- \_\_\_ 8. Compare the result of the scenario to that of a known good hearing aid

### Summary of Data

Were we able to meet the maximum amplification? Yes \_\_\_\_\_ No \_\_\_\_\_

Are modifications to the algorithm needed? Yes \_\_\_\_\_ No \_\_\_\_\_

Place Result Graph here

Testing Part 3 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 4 – Test for S4

- \_\_\_ 1. Flash the latest code to the evaluation board
- \_\_\_ 2. connect all of the subsystems together
- \_\_\_ 3. Ensure batteries are charged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 4. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 5. Place the microphones and speaker inside of the audiologist test equipment
- \_\_\_ 6. Run a scenario that will test at the limits of the frequency ranges to see hearing aid response
- \_\_\_ 7. Compare the result of the scenario to that of a known good hearing aid

### Summary of Data

Were we able to amplify the necessary frequencies? Yes \_\_\_\_\_ No \_\_\_\_\_

Are modifications to the algorithm needed? Yes \_\_\_\_\_ No \_\_\_\_\_

Place Result Graph here

Testing Part 4 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 5 – Test for S7

- \_\_\_ 1. Connect all of the subsystems together
- \_\_\_ 2. Ensure batteries are uncharged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 3. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 4. Plug the USB Cable into the battery charger
- \_\_\_ 5. Measure the battery voltage as the charging begins

### Summary of Data

Did the Batteries start charging? Yes \_\_\_\_\_ No \_\_\_\_\_

Did they hold a charge Yes \_\_\_\_\_ No \_\_\_\_\_

Testing Part 5 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 6 – Test for S5

- \_\_\_ 1. connect all of the subsystems together
- \_\_\_ 2. Ensure batteries are uncharged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 3. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 4. Plug the USB Cable into the battery charger
- \_\_\_ 5. Measure the battery voltage as the charging begins, start timing
- \_\_\_ 6. When batteries are fully charged stop timing. Record time to charge

### Summary of Data

Was the time to charge acceptable? Yes \_\_\_\_\_ No \_\_\_\_\_

What was the time to charge? \_\_\_\_\_

Testing Part 6 Sign Off \_\_\_\_\_ Date \_\_\_\_\_

# Hearing Aid Redesign: Test Plans

## EE 7: System Test

Date Completed \_\_\_\_\_

Performed By \_\_\_\_\_

### Part 7 – Test for S12

- \_\_\_ 1. Connect all of the subsystems together
- \_\_\_ 2. Ensure batteries are fully charged, charger circuit is correct, and amplification and filtering is in place
- \_\_\_ 3. Check that the microphones and speaker are connected, as well as the switch.
- \_\_\_ 4. Set hearing aid to a mode that will use the most power
- \_\_\_ 5. Measure the battery voltage and plug them into the hearing aid
- \_\_\_ 6. Start timing, and monitor the battery voltage
- \_\_\_ 7. When the batteries are fully drained stop timing
- \_\_\_ 8. Record the worst case battery life

### Summary of Data

Was the battery life of the hearing aid acceptable? Yes \_\_\_\_\_ No \_\_\_\_\_

What was the battery life in this worst case? \_\_\_\_\_

Do changes need to be made to extend the battery life? Yes \_\_\_\_\_ No \_\_\_\_\_

If so list any possible changes here:

Testing Part 6 Sign Off \_\_\_\_\_ Date \_\_\_\_\_